



## Vacuum Desiccator Cabinet



**Price:** RFQ  
**Product Code:** 1300-1-A  
**Manufacturer:** CLEATECH LLC  
**Made in:** U.S.A

**Lead-time** Usually ship in 2-3 week  
**Description:** 1300 Series economical Vacuum Desiccator cabinets designed for vacuum testing, degassing small parts, defoaming, mixed liquids, biological tests and storing samples under sustained low and intermediate vacuums. These vacuum desiccators are available in standard and custom sizes.

### Standard Features:

- Vacuum Internal Dimensions 7"x7"x8"
- Made of clear 1" thickness clear acrylic material
- Removable cover ( Type-I ) makes it easy to insert and take out samples
- Vacuum gasket on door ensures airtight seal
- Supports vacuum level of 29.5" Hg and more ( > 756 millimeter of Mercury)
- Include one removable perforated Acrylic shelf
- Two vacuum valves mounted to top allow purging of the desiccator with nitrogen before vacuum is applied, 3/8" ID barbed fittings included
- Stainless steel latches completely corrosion resistance

### Standard Options:

- Vacuum Gauge 0 to 30 in. Hg ( Installed )
- Additional Perforated Acrylic Shelf
- Two-stage Continuous Duty, Rotary Vacuum Pump 5.5 CFM
- Tubing, Vacuum-Rated 3-A Sanitary Clear PVC, 3/8" ID, feet length
- Vacuum Control Unit, 120V Monitors the vacuum level and automatically maintains a preset vacuum level ( to 29.8" Hg )

### Package Info:

- Package Type: Crate Cardboard box
- Package Dimensions: 17"x17"x19"
- Approximate Gross Weight: 25 lbs.