



Vacuum Desiccator Cabinet



Price: RFQ
Product Code: 1300-2-C
Manufacturer: CLEATECH LLC
Made in: U.S.A

Lead-time Usually ship in 2-3 week
Description: Cleatech's economical Vacuum Desiccator cabinets designed for vacuum testing, degassing small parts, defoaming, mixed liquids, biological tests and storing samples under sustained low and intermediate vacuums. These vacuum desiccators are available in standard and custom sizes.

Standard Features:

- Vacuum Internal Dimensions 10"x12"x16"
- Made of clear 1" thickness clear acrylic material
- Front Hinged Swing Door (Type-II) for easier access to the samples in all shelves
- Vacuum gasket on door ensures airtight seal
- Supports vacuum level of 29.8" Hg and more (> 756 millimeter of Mercury)
- Include one removable perforated Acrylic shelf
- Two vacuum valves mounted to top allow purging of the desiccator with nitrogen
Before vacuum is applied, 3/8" ID barbed fitting included
- Stainless steel latches completely corrosion resistance

Standard Options:

- Vacuum Gauge 0 to 30 in. Hg (Installed)
- Additional Perforated Acrylic Shelf
- Two-stage Continuous Duty, Rotary Vacuum Pump 5.5 CFM
- Tubing, Vacuum-Rated 3-A Sanitary Clear PVC, 3/8" ID, feet length
- Vacuum Control Unit, 120V Monitors the vacuum level and automatically maintains a preset vacuum level (to 29.8" Hg)

Package Info:

- Package Type: Crate
- Package Dimensions: 20"x20"x24"
- Approximate Gross Weight: 60 lbs.