



## **Stainless Steel Desiccator – Dry storage Cabinet- 1600 Series**



**Price:** RFQ

**Product Code:** 1600-1-H

**Manufacturer:** CLEATECH LLC

**Made in:** U.S.A

**Lead-time** Usually ship in 2-3 week

**Description:** Stainless steel desiccators (dry boxes) are durable, low humidity storage cabinets which designed for maximum weight capacity and long term durability.

Each Stainless Steel desiccator is a Non-vacuum dry cabinet can be equipped with optional inlet and outlet ports that allow gas to be purged through the storage area. For automated Nitrogen purging and humidity control you can use Cleatech's Automatic purge control unit. Custom sizes and features are available as required.

### **Standard Features:**

- Stainless Steel Desiccator Overall Dimensions: 22"Wx 22"Dx 16"H
- Usable Shelf Space: 21 ½"W x 21 ½"D
- Bench-top, Single Chamber Desiccator cabinet
- 304 Stainless steel Gage 16 heavy-duty construction
- Stainless Steel door hinge for maximum door support
- One Chrome Plated Lift-and-turn Compression Latch for perfect door Seals
- One-piece non-adhesive door sealing
- Chrome-plated rack allow shelf placement at one-inch increments
- 304 stainless steel door frame with Clear static-dissipative PVC window

### **Standard Options:**

- One-Channel Automatic Nitro Purge Control Unit (PCU) with one humidity sensor.
- Flowmeter, Range 2-20 SCFH With Stainless steel Valve
- Gas-In Port & Adjustable Auto pressure relief / Bleeding valve for Nitrogen / Argon purge
- Heavy duty 304 Stainless steel perforated shelves for upto 45 lbs. weight
- One to Four level Heavy-duty Power coated steel stand
- Electropolished Exterior & Passivated interior for cleanroom compatibility
- Cam lock, Same Key, installed for each door / Cam lock, Different Key.

### **Package Info:**

- Package Type: Crate
- Package Dimensions: 28"x28"x22"
- Approximate Gross Weight: 66 lbs.